


## Product Datasheet

### CCP110 antibody (orb40845)

|                            |                                                                               |
|----------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>         | CCP110 antibody                                                               |
| <b>Species/Host</b>        | Rabbit                                                                        |
| <b>Reactivity</b>          | Human                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                  |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                                |
| <b>Immunogen</b>           | Recombinant Human Centriolar coiled-coil protein of 110 kDa protein (1-280AA) |
| <b>Target</b>              | CCP110                                                                        |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                             |
| <b>Form/Appearance</b>     | Liquid                                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                 |
| <b>Note</b>                | For research use only                                                         |
| <b>Isotype</b>             | IgG                                                                           |
| <b>Clonality</b>           | Polyclonal                                                                    |
| <b>Uniprot ID</b>          | <b>O43303</b>                                                                 |
| <b>Dilution Range</b>      | IHC-P:1:20-200                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                               |



Western blot analysis of HeLa(lane 1), 2...