

Product Datasheet

Mouse OSCAR ELISA Kit (orb408451)

Description	Mouse Osteoclast-associated immunoglobulin-like receptor(OSCAR) ELISA kit
Reactivity	Mouse
Range	37.5 pg/mL-2400 pg/mL
Target	OSCAR
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8VBT3
Sensitivity	9.375 pg/mL
Expiration Date	Please enquire.