

Product Datasheet

Mouse PAH ELISA Kit (orb408466)

Description	Mouse Phenylalanine-4-hydroxylase(PAH) ELISA kit
Reactivity	Mouse
Range	23.5 pg/mL-1500 pg/mL
Target	PAH
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P16331
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.