

Product Datasheet

Mouse GPX4 ELISA Kit (orb408474)

Description	Mouse Phospholipid hydroperoxide glutathione peroxidase, mitochondrial(GPX4) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	GPX4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q91XR9
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.