

Product Datasheet

Mouse SERPING1 ELISA Kit (orb408478)

Description	Mouse Plasma protease C1 inhibitor(SERPING1) ELISA kit
Reactivity	Mouse
Range	5.6 ng/mL-360 ng/mL
Target	SERPING1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P97290
Sensitivity	1.4 ng/mL
Expiration Date	Please enquire.