

Product Datasheet

Mouse PDGF-BB ELISA Kit (orb408486)

Description	Mouse Platelet-derived growth factor subunit BB(PDGF-BB) ELISA kit
Reactivity	Mouse
Range	93.75 pg/mL-6000 pg/mL
Target	PDGF-BB
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Sensitivity	23.44 pg/mL
Expiration Date	Please enquire.