

Product Datasheet

Mouse PRDM16 ELISA Kit (orb408491)

Description	Mouse PR domain zinc finger protein 16(PRDM16) ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	PRDM16
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	A2A935
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.