

Product Datasheet

Mouse PTGDS ELISA Kit (orb408505)

Description	Mouse Prostaglandin-H2 D-isomerase(PTGDS) ELISA kit
Reactivity	Mouse
Range	4.7 ng/mL-300 ng/mL
Target	PTGDS
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	009114
Sensitivity	1.2 ng/mL
Expiration Date	Please enquire.