

Product Datasheet

Mouse PARK7 ELISA Kit (orb408518)

Description	Mouse Protein DJ-1(PARK7) ELISA kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	PARK7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q99LX0
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.