

Product Datasheet

Mouse NOTUM ELISA Kit (orb408528)

Description	Mouse Protein notum homolog(NOTUM) ELISA kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	NOTUM
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Uniprot ID	Q8R116
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.