

Product Datasheet

Mouse PADI4 ELISA Kit (orb408554)

Description	Mouse Protein-arginine deiminase type-4(PADI4) ELISA kit
Reactivity	Mouse
Range	4.7 pg/mL-300 pg/mL
Target	PADI4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9Z183
Sensitivity	1.17 pg/mL
Expiration Date	Please enquire.