

Product Datasheet

Mouse WNT3 ELISA Kit (orb408561)

Description	Mouse Proto-oncogene Wnt-3(WNT3) ELISA kit
Reactivity	Mouse
Range	39 pg/ml-2500 pg/ml
Target	WNT3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P17553
Sensitivity	9.75 pg/ml
Expiration Date	Please enquire.