

Product Datasheet

Mouse RIPK3 ELISA Kit (orb408571)

Description	Mouse Receptor-interacting serine/threonine-protein kinase 3(RIPK3) ELISA kit
Reactivity	Mouse
Range	46.88 pg/mL-3000 pg/mL
Target	RIPK3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9QZL0
Sensitivity	11.72 pg/mL
Expiration Date	Please enquire.