

Product Datasheet

Mouse CD163 ELISA Kit (orb408595)

Description	Mouse Scavenger receptor cysteine-rich type 1 protein M130(CD163) ELISA kit
Reactivity	Mouse
Range	0.93 ng/ml-60ng/ml
Target	CD163
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q2VLH6
Sensitivity	0.23 ng/ml
Expiration Date	Please enquire.