

Product Datasheet

Mouse SFRP4 ELISA Kit (orb408601)

Description	Mouse Secreted frizzled-related protein 4(SFRP4) ELISA kit
Reactivity	Mouse
Range	0.47 ng/mL-30 ng/mL
Target	SFRP4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9Z1N6
Sensitivity	0.11 ng/mL
Expiration Date	Please enquire.