

Product Datasheet

Mouse SCGB3A1 ELISA Kit (orb408607)

Description	Mouse Secretoglobin family 3A member 1(SCGB3A1) ELISA kit
Reactivity	Mouse
Range	62.5 pg/ml -4000 pg/ml
Target	SCGB3A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q920D7
Sensitivity	15.6 pg/ml
Expiration Date	Please enquire.