

Product Datasheet

Mouse SCGB3A2 ELISA Kit (orb408608)

Description	Mouse Secretoglobin family 3A member 2(SCGB3A2) ELISA kit
Reactivity	Mouse
Range	9.4 pg/mL-600 pg/mL
Target	SCGB3A2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q920H1
Sensitivity	2.35 pg/mL
Expiration Date	Please enquire.