

Product Datasheet

Mouse SEMA3D ELISA Kit (orb408616)

Description	Mouse Semaphorin-3D(SEMA3D) ELISA kit
Reactivity	Mouse
Range	23.5 pg/mL-1500 pg/mL
Target	SEMA3D
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8BH34
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.