

Product Datasheet

Mouse SCUBE1 ELISA Kit (orb408642)

Description	Mouse Signal peptide, CUB and EGF-like domain-containing protein 1(SCUBE1) ELISA kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	SCUBE1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q6NZL8
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.