

Product Datasheet

Mouse SOD2 ELISA Kit (orb408663)

Description	Mouse Superoxide dismutase [Mn], mitochondrial(SOD2) ELISA kit
Reactivity	Mouse
Range	0.625 ng/mL-40 ng/mL
Target	SOD2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P09671
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.