

Product Datasheet

Mouse TXNIP ELISA Kit (orb408689)

Description	Mouse Thioredoxin-interacting protein(TXNIP) ELISA kit
Reactivity	Mouse
Range	0.45 ng/ml-30 ng/ml.
Target	TXNIP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8BG60
Sensitivity	0.11 ng/ml
Expiration Date	Please enquire.