

Product Datasheet

Mouse Thrombospondin 2 ELISA Kit (orb408690)

Description	Mouse Thrombospondin-2 (THBS2) ELISA kit
Reactivity	Mouse
Range	6.25 ng/ml-400 ng/ml.
Target	THBS2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates.
Uniprot ID	Q03350
Sensitivity	1.56 ng/ml.
Expiration Date	Please enquire.