

Product Datasheet

Mouse TSLP ELISA Kit (orb408694)

Description	Mouse Thymic stromal lymphopoietin(TSLP) ELISA kit
Reactivity	Mouse
Range	15.6 pg/ml-1000 pg/ml
Target	TSLP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9JIE6
Sensitivity	3.9 pg/ml
Expiration Date	Please enquire.