



## Product Datasheet

### GNA15 antibody (orb40870)

|                            |                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------|
| <b>Description</b>         | GNA15 antibody                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                  |
| <b>Reactivity</b>          | Human                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                            |
| <b>Tested Applications</b> | ELISA, IHC                                                                              |
| <b>Immunogen</b>           | Recombinant Human Guanine nucleotide-binding protein subunit alpha-15 protein (1-374AA) |
| <b>Target</b>              | GNA15                                                                                   |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                                       |
| <b>Form/Appearance</b>     | Liquid                                                                                  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                           |
| <b>Note</b>                | For research use only                                                                   |
| <b>Isotype</b>             | IgG                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                              |
| <b>Uniprot ID</b>          | <b>P30679</b>                                                                           |
| <b>Dilution Range</b>      | IHC-P:1:20-200                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                         |



Immunohistochemical staining of human sk...



Immunohistochemical staining of human sm...