

Product Datasheet

Mouse TGFB β 3 ELISA Kit (orb408705)

Description	Mouse Transforming growth factor beta receptor type 3(TGFB β 3) ELISA kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	TGFB β 3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O88393
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.