

Product Datasheet

Mouse TPT1 ELISA Kit (orb408707)

Description	Mouse Translationally-controlled tumor protein(TPT1) ELISA kit
Reactivity	Mouse
Range	9.37 pg/mL-600 pg/mL
Target	TPT1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P63028
Sensitivity	2.34 pg/mL
Expiration Date	Please enquire.