

Product Datasheet

Mouse TNFSF13 ELISA Kit (orb408725)

Description	Mouse Tumor necrosis factor ligand superfamily member 13(TNFSF13) ELISA kit
Reactivity	Mouse
Range	93.75 pg/mL-6000 pg/mL
Target	TNFSF13
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	Q9D777
Sensitivity	23.43 pg/mL
Expiration Date	Please enquire.