

Product Datasheet

Mouse TNFSF4 ELISA Kit (orb408726)

Description	Mouse Tumor necrosis factor ligand superfamily member 4(TNFSF4) ELISA kit
Reactivity	Mouse
Range	15.6pg/mL-1000pg/mL
Target	TNFSF4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P43488
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.