

Product Datasheet

Mouse TNFSF9 ELISA Kit (orb408727)

Description	Mouse Tumor necrosis factor ligand superfamily member 9(TNFSF9) ELISA kit
Reactivity	Mouse
Range	23.5 pg/mL-1500 pg/mL
Target	TNFSF9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P41274
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.