

Product Datasheet

Mouse TNFRSF13C ELISA Kit (orb408731)

Description	Mouse Tumor necrosis factor receptor superfamily member 13C(TNFRSF13C) ELISA kit
Reactivity	Mouse
Range	39 pg/ml - 2500 pg/ml
Target	TNFRSF13C
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9D8D0
Sensitivity	9.75 pg/ml
Expiration Date	Please enquire.