

Product Datasheet

Mouse TNFRSF18 ELISA Kit (orb408734)

Description	Mouse Tumor necrosis factor receptor superfamily member 18(TNFRSF18) ELISA kit
Reactivity	Mouse
Range	78 pg/ml - 5000 pg/ml
Target	TNFRSF18
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O35714
Sensitivity	19.5 pg/ml
Expiration Date	Please enquire.