

Product Datasheet

Mouse WFDC2 ELISA Kit (orb408768)

Description	Mouse WAP four-disulfide core domain protein 2(WFDC2) ELISA kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	WFDC2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9DAU7
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.