

Product Datasheet

Mouse IP3 ELISA Kit (orb408774)

Description	Mouse inositol 1,4,5,-trisphosphate(IP3) ELISA Kit
Reactivity	Mouse
Range	0.78 pg/mL-50 pg/mL
Target	IP3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Sensitivity	0.195 pg/mL
Expiration Date	Please enquire.