

## Product Datasheet

### Mouse NTX ELISA Kit (orb408784)

|                          |                                                                   |
|--------------------------|-------------------------------------------------------------------|
| <b>Description</b>       | Mouse cross linked N-telopeptide of type I collagen,NTX ELISA Kit |
| <b>Reactivity</b>        | Mouse                                                             |
| <b>Range</b>             | 3.12 nM BCE-200 nM BCE                                            |
| <b>Target</b>            | NTX                                                               |
| <b>Note</b>              | For research use only                                             |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                                 |
| <b>Sensitivity</b>       | 0.78 nM BCE                                                       |
| <b>Expiration Date</b>   | Please enquire.                                                   |