

Product Datasheet

Mouse ADAMTS13 ELISA Kit (orb408795)

Description	Mouse ADAM metallopeptidase with thrombospondin type 1 motif, 13 (ADAMTS13) ELISA kit
Reactivity	Mouse
Range	78 pg/mL-5000 pg/mL
Target	ADAMTS13
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q769J6
Sensitivity	19.5 pg/mL
Expiration Date	Please enquire.