

## Product Datasheet

### Mouse CCL25 ELISA Kit (orb408804)

|                          |                                                  |
|--------------------------|--------------------------------------------------|
| <b>Description</b>       | Mouse C-C motif chemokine<br>25(CCL25) ELISA kit |
| <b>Reactivity</b>        | Mouse                                            |
| <b>Range</b>             | 78 pg/ml - 5000 pg/ml                            |
| <b>Target</b>            | CCL25                                            |
| <b>Note</b>              | For research use only                            |
| <b>Application notes</b> | Detection Wavelength: 450 nm                     |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                |
| <b>Uniprot ID</b>        | <b>O35903</b>                                    |
| <b>Sensitivity</b>       | 19.5 pg/ml                                       |
| <b>Expiration Date</b>   | Please enquire.                                  |