

Product Datasheet

Mouse Interferon gamma ELISA Kit (orb408817)

Description	Mouse C-X-C motif chemokine 11(CXCL11) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	CXCL11
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates,tissue homogenates
Uniprot ID	Q9JHH5
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.