

Product Datasheet

Mouse CXCL1 ELISA Kit (orb408819)

Description	Mouse chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha) (CXCL1) ELISA k...
Reactivity	Mouse
Range	7.8 pg/mL-500 pg/mL
Target	CXCL1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P12850
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.