

Product Datasheet

Mouse CXCL16 ELISA Kit (orb408820)

Description	Mouse CXC-chemokine ligand 16,CXCL16 ELISA Kit
Reactivity	Mouse
Range	7.8125 pg/mL-500 pg/mL
Target	CXCL16
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma,tissue homogenates
Uniprot ID	Q8BSU2
Sensitivity	3.012 pg/mL
Expiration Date	Please enquire.