

Product Datasheet

Mouse LDHD ELISA Kit (orb408828)

Description	Mouse D-Lactate Dehydrogenase,D-LDH ELISA Kit
Reactivity	Mouse
Range	4.69 mU/mL-300 mU/mL
Target	LDHD
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q7TNG8
Sensitivity	1.17 mU/mL
Expiration Date	Please enquire.