

Product Datasheet

Mouse IL4I1 ELISA Kit (orb408840)

Description	Mouse L-amino-acid oxidase(IL4I1) ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	IL4I1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O09046
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.