

Product Datasheet

Mouse SIRT1 ELISA Kit (orb408848)

Description	Mouse NAD-dependent deacetylase sirtuin-1(Sirt1)ELISA Kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	SIRT1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q923E4
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.