

Product Datasheet

Mouse NAG ELISA Kit (orb408851)

Description	Mouse N-acetyl- β -D-glucosaminidase,NAG ELISA Kit
Reactivity	Mouse
Range	15.6 μ IU/mL-1000 μ IU/mL
Target	NAG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P20060
Sensitivity	3.9 μ IU/mL
Expiration Date	Please enquire.