

Product Datasheet

Mouse WNT7A ELISA Kit (orb408877)

Description	Mouse Protein Wnt-7a(WNT7A) ELISA kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	WNT7A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P24383
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.