

Product Datasheet

Mouse B2M ELISA Kit (orb408892)

Description	Mouse β 2-microglobulin,BMG/ β 2-MG ELISA Kit
Reactivity	Mouse
Range	0.156 μ g/ml-10 μ g/ml.
Target	B2M
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, urine, tissue homogenates
Uniprot ID	P01887
Sensitivity	0.156 μ g/ml.
Expiration Date	Please enquire.