

## Product Datasheet

### Mouse HEXA ELISA Kit (orb408894)

|                          |                                                              |
|--------------------------|--------------------------------------------------------------|
| <b>Description</b>       | Mouse $\beta$ -hexosaminidase A, $\beta$ -Hex A<br>ELISA Kit |
| <b>Reactivity</b>        | Mouse                                                        |
| <b>Range</b>             | 3.12 $\mu$ U/mL-200 $\mu$ U/mL                               |
| <b>Target</b>            | HEXA                                                         |
| <b>Note</b>              | For research use only                                        |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                 |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates                         |
| <b>Uniprot ID</b>        | <b>P29416</b>                                                |
| <b>Sensitivity</b>       | 0.78 $\mu$ U/mL                                              |
| <b>Expiration Date</b>   | Please enquire.                                              |