

Product Datasheet

Mouse MSH Alpha ELISA Kit (orb408908)

Description	Mouse Pro-opiomelanocortin (POMC) ELISA KIT
Reactivity	Mouse
Range	0.73 pg/mL-3000 pg/mL
Target	POMC
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P01193
Sensitivity	0.18 pg/mL
Expiration Date	Please enquire.