

Product Datasheet

Mouse CEACAM5 ELISA Kit (orb408911)

Description	Mouse carcinoembryonic antigen,CEA ELISA Kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	CEACAM5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q3UKK2
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.