

Product Datasheet

Mouse IL1RA ELISA Kit (orb408933)

Description	Mouse Interleukin 1 receptor antagonist, IL-1ra ELISA kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	IL1RN
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P25085
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.