

Product Datasheet

Mouse C19orf10 ELISA Kit (orb408939)

Description	Mouse interleukin 25(IL25/C19orf10) ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	C19orf10
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9CPT4
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.